



NEW

NF 3000R

HIGH CAPACITY REFRIGERATED BENCH TOP CENTRIFUGE

The NF 3000R high capacity refrigerated benchtop centrifuge has been designed with the needs of large volume laboratories and blood banks in mind.

Versatile Use

- Possibility to perform many centrifuge operations with different rotors
- Opening rotor with 4x750 ml capacity, offering different gode options for the most used tubes.

Opportunity to centrifuge 4 triple blood bags

Product Protection

- Programmable temperature between -10°C / +40°C and powerful cooling system that provides +4°C inside the cell at maximum speed for the safety of biologically active samples

Smart Control System

- Intelligent control system with 7" color touch screen
- User-friendly control system with easy-to-understand icons
- 3GB internal memory for centrifugation processes
- USB port for recording of running programs
- Ability to define 3 users (Admin, Operator1, Operator2)
- Authorizing users

WIRELESS COMMUNICATION

- * N-Mobile™ and NuveControl™ remote control
- * Work records of each user
- * Warning system in case of any error



Max Speed	Swing Out Rotor: 4,600 rpm Angle Rotor: 15,000 pm
Maximum RCF	Swing Out Rotor: 4.637xg Angle Rotor: 28.677xg
Tube Capacity	Swing Out Rotor:4x750 ml / 4xTriple Blood Bag Angle Rotor: 12x100 ml
Control System	Angle Rotor: 30x15 ml Angle Rotor: 15x10 ml Angle Rotor: 72x1.5/2 ml
Indicator	N-Wise™ Programmable Microprocessor Control System
Speed Adjustment Range	7" Color Touch Screen
Speed Adjustment Range	500-15,000rpm
Speed Control Accuracy	10rpm
Timer Setting Range	± 20 rpm
Timer Setting Range	1-99 minutes and indefinitely
Acceleration/Braking Range	1 min
Temperature Setting Range	10/10
Temperature Adjustment Level	-10°C / +40°C
Engine	1°C
Body and Cover Structure	Maintenance-free induction motor
Chamber Structure	Epoxy-polyester powder coated steel
Installed Power	304 stainless steel



NÜVE SANAYİ MALZEMELERİ İMALAT VE TİCARET A.Ş.

Saracalar Mah. Saracalar Kümeevleri No: 4/2 Akyurt 06750 ANKARA / TÜRKİYE t. +90 312 399 28 30 f. +90 312 399 21 97

nuve.com.tr satis@nuve.com.tr

ISO 9001: 2015
EN ISO 13485: 2016

